


How Do I Use a Tablet Splitter?


Not all medications are appropriate for splitting. Talk to your pharmacist if you have questions.

Step-by-Step Instructions:

1. Place tablet splitter on a flat, stable surface.
2. Open the tablet splitter. Notice the blade attached on the top as you can see in the picture below. DO NOT TOUCH the blade. It is very sharp.



3. Place tablet within the V-shaped holder. As shown in the diagram below, try to place the tablet so the blade will cut it into two equal pieces. If there is a score in the middle of the tablet, line it up with the blade.



4. Using one smooth and steady motion, close the top down to allow the tablet to split.
5. Open and remove the tablet halves. We suggest you split one day's worth at a time instead of splitting the entire supply at once. You may temporarily store unused portion of the day's supply within storage compartments of the splitter.
6. Always store your tablet splitter closed and out of reach of children.



- Cleaning Instruction: Wipe off any leftover powder with a tissue as necessary. You may rinse the splitter with warm soapy water and allow it to air dry before your next use. Do not wash in dishwasher.
- Physicians Plus will provide one tablet splitter for each member of household upon request.

P+5322-0904